

Accusol

Dialysis solution for CRRT



Clear and simple

ACCUSOL dialysis solution for CRRT offers healthcare professionals a 100 % bicarbonate buffered solution which is completely lactate-free.

Available in a range of potassium chloride concentrations: 0, 2 or 4 mmol/l

Environmentally friendly: PVC, DEHP and Latex free presentation

Ease of use

- Easy, two-step mixing process
- Double seal ensures fluid cannot be infused unless mixed
- Simple Luer connection – no frangible seal
- Colour-coding to aid potassium strength recognition


One buffer, 100 % bicarbonate, lactate free

Two-step quick mixing

Three colour-coded formulations

ACCUSOL

Dialysis solution for CRRT



	ACCUSOL K ⁺ 0 (Potassium 0 mmol/l)	ACCUSOL K ⁺ 2 (Potassium 2 mmol/l)	ACCUSOL K ⁺ 4 (Potassium 4 mmol/l)
COMPOSITION (mmol/l)*			
Upper chamber			
Sodium	129	129	129
Potassium	0	2.67	5.33
Calcium	2.33	2.33	2.33
Magnesium	0.667	0.667	0.667
Chloride	145.7	148.4	151.1
Glucose	0	7.4	7.4
Lower chamber			
Sodium	173	173	173
Bicarbonate	160	160	160
Ready to use solution			
Sodium	140	140	140
Potassium	0	2	4
Calcium	1.75	1.75	1.75
Magnesium	0.5	0.5	0.5
Chloride	109.3	111.3	113.3
Bicarbonate	35	35	35
Glucose	0	5.55	5.55
Osmolarity (mOsm/L)	287	296	300
DESCRIPTION			
Product Code	RPE9701	RPE9702	RPE9703
Volume	5000 ml	5000 ml	5000 ml
Units per case	2	2	2

* Converted using formula mol = mass / Relative Atomic Mass
 Figures based on Direction Insert

- **Individualised prescribing.** Three ready-to-mix formulations offering multiple potassium concentrations, better matching your patient prescription
- **100% bicarbonate buffer.** ACCUSOL dialysis solution contains the necessary electrolytes to ensure physiologic pH with zero lactate
- **Integral Medication port.** Enabling potassium chloride supplementation as the patient requires
- **Colour-coded.** Potassium concentrations are identified by a colour-coded hanger hole
- **Confidence in solution flow.** The double seal and the bag connection system with wings ensure unimpeded flow and reduce the risk of infusing unmixed solution

CRRT (Continuous Renal Replacement Therapy), PVC (Polyvinyl chloride), DEHP (Diethylhexyl phthalate).

ACCUSOL Solution*
 Nikkiso Belgium bvba
 Industriepark 6
 3300 Tienen
 Belgium
 CE 0123

www.nikkisoabp.com